



Packing HiScale™, XK, and Tricorn™ chromatography columns with Capto™ and MabSelect™ resins

Table of contents

Impact of column packing on results and efficiency	3
Terminology used in this guide	4
Slurry concentration	4
Compression	4
Custom Design Media (CDM) group	4
Planning to scale up?	4
Column packing recommendations	5
Slurry preparation	5
Use a packing tube	5
Test the efficiency of the packed column	5
Not enough time to pack your own columns?	6
Packing parameters by chromatography resin	7
Ion exchange chromatography	
Capto Q ImpRes, Capto SP ImpRes, Capto S ImpAct	7
Capto Q, Capto S, Capto DEAE	8
Affinity chromatography	
Capto Heparin	8
MabSelect PrismA, MabSelect	9
MabSelect Xtra, MabSelect SuRe™, MabSelect SuRe LX, Capto L	10
VIII Select	10
KappaSelect, Lambda FabSelect, IgG Select, Capto Blue, Capto Blue (High Sub), Capto Chelating, VII Select, CaptoDeViris, and Alpha-1-Antitrypsin Select	14
Multimodal chromatography	
Capto Core 400	8
Capto Core 700, Capto MMC	11
Capto adhere, Capto adhere ImpRes	12
Capto MMC ImpRes	13
Hydrophobic interaction chromatography	
Capto Phenyl ImpRes, Capto Butyl ImpRes	13
Capto Phenyl (High Sub), Capto Butyl, Capto Octyl	14

Introduction

To obtain optimal performance during your protein purification, it is critical to use well-packed columns. For an optimal column packing, there are several parameters that need to be controlled during packing. These packing parameters are specific to the resin, and the type and dimension of the column that needs to be packed.

This guide summarizes parameters for packing small-scale HiScale, XK, and Tricorn columns with Capto and MabSelect resins from GE Healthcare Life Sciences.

Impact of column packing on results and efficiency

The way that liquid flows through a chromatography column depends on how it is packed. A well-packed bed generates a stable column that offers good resolution. An even flow of liquid through the column will generate peaks that are narrow and sharp, as shown in Fig 1B.




Fig 1. Consequence of a well packed resin bed: (A) Even flow of liquid; (B) Narrow and sharp peaks on the chromatogram.

A poorly packed bed generates an uneven flow through the column. Peaks will be broad, as shown in Fig 2B.




Fig 2. Consequence of a poorly packed resin bed: (A) Uneven flow of liquid; (B) Broad peaks on the chromatogram.

A well-packed column ensures reproducible column performance run after run. When liquid moves unevenly through the column, protein peaks will be broad and might overlap with other neighboring peaks. Therefore, the purity of the target protein could be negatively impacted when using a poorly packed column.

Terminology used in the guide

Slurry concentration

The slurry concentration is used to facilitate the calculation of the amount of resin needed to pack a certain bed height. For successful packing, different resins require different slurry concentrations. To determine the slurry concentration, we recommend that you use the [Slurry Concentration Kit](#). You can also watch our video [Column packing tutorial: Determining slurry concentration using a slurry kit](#).

Other methods involving centrifugation and sedimentation are also possible to use.

Note: The level of accuracy for determining the [slurry concentration](#) is not as critical as it is when packing large-scale columns. Small-scale columns may be packed with an excess of resin, which is removed after packing.

Compression

Depending on the column type and resin, the compression may be performed in two different ways:

- Mechanical compression with packing factor: after applying the settling flow, note the height of the consolidated bed before stopping the flow (after the flow is stopped the bed can slightly expand). The final bed height is calculated by dividing the consolidated bed height with the packing factor (PF):

Final bed height = Consolidated bed height/packing factor (PF)

Set the adapter against the consolidated packed bed, tighten the O-ring, and turn the end cap down until the calculated final bed height is reached.

- Flow compression: after applying the packing flow, the adapter is moved a specific distance (mm) into the packed bed to avoid gap formation.

Custom Design Media (CDM) group

If the chromatography resin you require is not part of our standard product offering, our specialists can custom design your chromatography resin from an appropriate base matrix and ligand.

Please note that CDM products are not supplied with instruction manuals.

Planning to scale up?

If you are working in process development and intend to scale-up, consider from the start the bed heights and flow rates that you will require at larger scales.

For consistent results when scaling up, we recommend that you keep residence time and/or bed heights constant.

Also, please note these packing methods are adapted to small-scale columns (< 5 cm i.d.) and are not always scalable to large-scale columns. If you need to scale-up, refer to the [packing recommendations](#) for large-scale Axichrom™, BPG, and Chromaflow™ columns.

Column packing recommendations

Slurry preparation

Ensure that the resin is washed thoroughly into the packing solution before starting the packing. Use a glass filter to wash the resin over to the packing solution. Suspend the resin by shaking and pour into the funnel and wash according to the following instructions:

1. Wash 5 times with 2 column volumes (CV) of packing solution. Gently stir with a plastic spatula between additions.
2. Pour the washed resin from the funnel into a beaker.
3. Add packing solution to obtain the indicated slurry concentration for specific resin and column type.

Use a packing tube

When packing HiScale, XK, and Tricorn, we recommend you use a packing connector together with an extra glass tube to serve as a packing reservoir (Fig 3). By using the packing tube, you can add the whole slurry volume in a single pouring. For detailed description on how to pack each column, we recommend that you download the column instructions document for the respective column types:




Fig 3. Connection of the packing tube for HiScale, XK, and Tricorn columns.

Test the efficiency of the packed column

Column efficiency testing plays a central role in the qualification and monitoring of packed bed performance. The desirable high column efficiency gives low band/peak broadening and is an indicator of how well packed the column is before starting purification.

Column efficiency is typically defined in terms of two parameters (Fig 4):

- Peak broadening over the column is described by an equivalent number of theoretical plates.
- Peak symmetry is described by a peak asymmetry factor, A_s

The procedure is described in each GE resin protocol and is also described in this Column efficiency testing application note [28937207](#).




Fig 4. Determination of column efficiency by number of theoretical plates per meter (N/m) and peak asymmetry factor (A_s).

Not enough time to pack your own columns?

Packing a column can be both difficult and time-consuming and needs practice. For SEC, columns can take up to a full day to pack. For other techniques, packing can require over an hour (see figure).

Instead of spending hours of your time packing and testing your column before you can start to purify your protein, consider purchasing a manufactured prepacked column. This eliminates the need for repeated column packing and ensures that you have a well-packed column for immediate use.

If the prepacked column you need is not available as a regular product, contact our [Custom Products](#) group. We can customize columns and resins to your specific requirements.

	Packing your own columns		Buying a prepacked column
	SEC	HIC/IEX/AC	Various chromatography techniques
Prepare slurry	~ 1 h	5 min	-
Prepare column	15-30 min	15-30 min	-
Pack the column	1-4 h	~ 30 min	-
Test the column	45-120 min	30 min	-
Time until column is ready to use	3-7.5 h	60-90 min	Column ready to use!

Capto Q ImpRes, Capto SP ImpRes, and Capto S ImpAct

Packing parameters

Your parameters				Resin preparation		Column packing						Post-column packing	
Resin	Packed bed height (cm)	Column type	Bed volume (mL)	Slurry/packing solution	Slurry concentration (%)	Settling flow The flow that consolidates the resin suspension (slurry)		Packing flow The flow that compresses the bed		Compression Bed is further compressed by turning the adapter manually into the bed		Conditioning flow The flow to get a uniform bed over the entire column height	
						mL/min	cm/h	mL/min	cm/h	Mechanical compression packing factor (PF) ¹	Flow compression (mm)	mL/min	cm/h
Capto Q ImpRes Capto SP ImpRes Particle size d_{50}^2 : 40 μm	10	Tricorn 5	2.0	10 mM NaCl	45–55	7.4	2250	7.4	2250	N/A	1	N/A	N/A
		Tricorn 10	7.9			29.4		29.4					
		HiScale 10	7.9	20% ethanol with 0.4 M NaCl	63	10.3	783	14.5	1108	1.05	N/A	14.5	1108
		XK 16 and HiScale 16	20.1	10 mM NaCl	45–55	40	1200	N/A	N/A	1.12		40	1200
		XK 26 and HiScale 26	53.1			106				1.12		106	
		XK 50 and HiScale 50	196.3			393				1.12		393	
	HiScale 10	15.7	20% ethanol with 0.4 M NaCl			63				10.3		783	
	20	XK 16 and HiScale 16	40.2	10 mM NaCl	45–55	33	1000	N/A	N/A	1.12		33	1000
		XK 26 and HiScale 26	106.1			88				1.12	88		
		XK 50 and HiScale 50	392.5			327				1.12	327		
	25	HiScale 10	19.6	20% ethanol with 0.4 M NaCl	63	10.3	783	14.5	1108	1.04	14.5	1108	
	35	XK 16 and HiScale 16	70.3	10 mM NaCl	45–55	23	700	N/A	N/A	1.09	23	700	
		XK 26 and HiScale 26	185.7			62				1.09	62		
		XK 50 and HiScale 50	686.9			229				1.09	229		
	Capto S ImpAct Particle size d_{50}^2 : 50 μm	10	Tricorn 5	2.0	0.4 M NaCl	45–55	3.5	1070	3.5	1070	N/A	1	N/A
Tricorn 10			7.9	14			14						
XK 16 and HiScale 16			20.1	20% ethanol with 0.2 M sodium acetate	45–55	10	300	25	750	N/A	N/A	N/A	N/A
XK 26 and HiScale 26		53.1	27			53		600					
XK 50 and HiScale 50		196.3	98			196		600					
20		XK 16 and HiScale 16	40.2			7.4	220	16.7	500				
		XK 26 and HiScale 26	106.1			19		35	400				
		XK 50 and HiScale 50	392.5			72		98	300				
35		XK 16 and HiScale 16	70.3	3.3	100	8.4	250						
		XK 26 and HiScale 26	185.7	9		22	250						
		XK 50 and HiScale 50	686.9	33		65	200						

¹ The PF values in this table are specific for the lab-packing method used. The PF shown here bears no relation to the PF described for large-scale column packing procedures.

² d_{50}^2 = median particle size of the cumulative volume distribution.

Capto Q, Capto S, Capto DEAE, Capto Core 400, and Capto Heparin

Packing parameters

Your parameters				Resin preparation		Column packing						Post-column packing	
Resin	Packed bed height (cm)	Column type	Bed volume (mL)	Slurry/packing solution	Slurry concentration (%)	Settling flow The flow that consolidates the resin suspension (slurry)		Packing flow The flow that compresses the bed		Compression Bed is further compressed by turning the adapter manually into the bed		Conditioning flow The flow to get a uniform bed over the entire column height	
						mL/min	cm/h	mL/min	cm/h	Mechanical compression packing factor (PF) ²	Flow compression (mm)	mL/min	cm/h
Capto Q Capto S Particle size d_{50v}^3 : 90 μm	10	Tricorn 5	2.0	10 mM NaCl	40–60	1.8	540	9.8	3000	N/A	1	N/A	N/A
		Tricorn 10	7.9			7.1		39.3					
		XK 16 and HiScale 16	20.1			25		66					
	XK 26 and HiScale 26	53.1	66	1.15	66								
	XK 50 and HiScale 50	196.3	250	1.15	250								
	XK 16 and HiScale 16	40.2	25	1.1	25								
	XK 26 and HiScale 26	106.1	66	1.1	66								
	XK 50 and HiScale 50	392.5	250	1.15	250								
	20	XK 16 and HiScale 16	70.3	25	1.06	14	420						
		XK 26 and HiScale 26	185.7	66	1.06	37							
		XK 50 and HiScale 50	686.9	250	1.1	140							
	Capto DEAE Capto Heparin¹ Capto Core 400¹ Particle size d_{50v}^3 : 90 μm	10	Tricorn 5	2.0	20% ethanol with 0.2 M NaCl	45–60	5	1500	5	1500	N/A	1	N/A
Tricorn 10			7.9	20			20						
XK 16 and HiScale 16			20.1	25			66		1.1				
XK 26 and HiScale 26		53.1	66	1.15	66								
XK 50 and HiScale 50		196.3	250	1.15	250								
XK 16 and HiScale 16		40.2	25	1.1	25								
XK 26 and HiScale 26		106.1	66	1.1	66								
XK 50 and HiScale 50		392.5	250	1.15	250								
20		XK 16 and HiScale 16	70.3	25	1.06	14	420						
		XK 26 and HiScale 26	185.7	66	1.06	37							
		XK 50 and HiScale 50	686.9	250	1.1	140							

¹ Part of the CDM program from GE. Packing protocols for CDM-developed resins are not verified but are based on packing protocols for standard products having similar pressure-flow characteristics.

² The PF values in this table are specific for the lab-packing method used. The PF shown here bears no relation to the PF described for large-scale column packing procedures.

³ d_{50v} = median particle size of the cumulative volume distribution.

MabSelect Prisma and MabSelect

Packing parameters

Your parameters				Resin preparation		Column packing						Post-column packing						
Resin	Packed bed height (cm)	Column type	Bed volume (mL)	Slurry/packing solution	Slurry concentration (%)	Settling flow The flow that consolidates the resin suspension (slurry)		Packing flow The flow that compresses the bed		Compression Bed is further compressed by turning the adapter manually into the bed		Conditioning flow The flow to get a uniform bed over the entire column height						
						mL/min	cm/h	mL/min	cm/h	Mechanical compression packing factor (PF) ¹	Flow compression (mm)	mL/min	cm/h					
MabSelect Prisma Particle size d _{50v} ² : 60 µm	10	Tricorn 5	2.0	20% ethanol with 0.2 M NaCl	45-55	1.96	600	1.96	600	N/A	1	N/A	N/A					
		Tricorn 10	7.9			7.85		7.85										
		HiScale 10	7.9			2.62		2.62						200	1.1	N/A	9.16	700
		XK 16 and HiScale 16	20.1			6.7		N/A						N/A	1.12		16.8	500
		XK 26 and HiScale 26	53.1			17.7									1.15		44.2	
		XK 50 and HiScale 50	196.3			65.4									1.16		114.5	
	HiScale 10	15.7	2.62	2.62	200	1.1	5.24		400									
	20	XK 16 and HiScale 16	40.2	6.7	200	N/A	N/A	1.1	13.4	400								
		XK 26 and HiScale 26	106.1	17.7				1.1	35.4									
		XK 50 and HiScale 50	392.5	65.4				1.14	98.2	300								
		HiScale 10	19.6	2.62				2.62	200	1.1	4.19	320						
	35	XK 16 and HiScale 16	70.3	6.7	45-55	N/A	N/A	1.1	7.7	230								
		XK 26 and HiScale 26	185.7	17.7				1.1	20.4									
		XK 50 and HiScale 50	686.9	65.4				1.1	75.3									
Tricorn 5		2.0	2.29	700				2.29	700		N/A	1	N/A	N/A				
Tricorn 10	7.9	9.16	9.16															
MabSelect Particle size d _{50v} ² : 85 µm	10	XK 16 and HiScale 16	20.1	20% ethanol with 0.4 M NaCl	45-55	10	300	N/A	N/A	N/A	N/A	25	750					
		XK 26 and HiScale 26	53.1			27						1.15		26				
		XK 50 and HiScale 50	196.3			100						1.15		250				
		XK 16 and HiScale 16	40.2			10						1.1		15				
	20	XK 26 and HiScale 26	106.1			27						300	N/A	N/A	1.13	40	450	
		XK 50 and HiScale 50	392.5			100						1.1	150					
		XK 16 and HiScale 16	70.3			10						1.06	8.6					
	35	XK 26 and HiScale 26	185.7			27						1.1	23	260				
		XK 50 and HiScale 50	686.9			100						1.06	86					

¹ The PF values in this table are specific for the lab-packing method used. The PF shown here bears no relation to the PF described for large-scale column packing procedures.

² d_{50v} = median particle size of the cumulative volume distribution.

MabSelect Xtra, VIII Select, MabSelect SuRe, MabSelect SuRe LX, and Capto L

Packing parameters

Your parameters				Resin preparation		Column packing						Post-column packing				
Resin	Packed bed height (cm)	Column type	Bed volume (mL)	Slurry/packing solution	Slurry concentration (%)	Settling flow The flow that consolidates the resin suspension (slurry)		Packing flow The flow that compresses the bed		Compression Bed is further compressed by turning the adapter manually into the bed		Conditioning flow The flow to get a uniform bed over the entire column height				
						mL/min	cm/h	mL/min	cm/h	Mechanical compression packing factor (PF) ²	Flow compression (mm)	mL/min	cm/h			
MabSelect Xtra VIII Select¹ Particle size d_{50v}^3 : 75 μm	10	Tricorn 5	2.0	0.2 M NaCl	25–50	0.5	150	3	900	N/A	1	N/A	N/A			
		Tricorn 10	7.9			2		12								
		XK 16 and HiScale 16	20.1			10		27						1.1	25	750
	XK 26 and HiScale 26	53.1	100	1.15	66											
	20	XK 50 and HiScale 50	196.3	20% ethanol with 0.4 M NaCl	45–55	10	300	N/A	N/A	1.1	N/A	15	450			
		XK 16 and HiScale 16	40.2			27								1.13	40	
		XK 26 and HiScale 26	106.1			100								1.15	150	
	35	XK 50 and HiScale 50	392.5	20% ethanol with 0.4 M NaCl	45–55	10	300	N/A	N/A	1.06	N/A	8.6	260			
		XK 16 and HiScale 16	70.3			27								1.1	23	
		XK 26 and HiScale 26	185.7			100								1.06	86	
	MabSelect SuRe MabSelect SuRe LX Capto L Particle size d_{50v}^3 : 90 μm	10	Tricorn 5	2.0	20% ethanol with 0.2 M NaCl	45–55	2.29	700	2.29	700	N/A	1	N/A	N/A		
			Tricorn 10	7.9			9.16		9.16							
HiScale 10			7.9	16			1222		16						1222	1.0
20		XK 16 and HiScale 16	20.1	20% ethanol with 0.4 M NaCl	45–55	10	300	N/A	N/A	1.1	N/A	15	450			
		XK 26 and HiScale 26	53.1			27								1.15	66	
		XK 50 and HiScale 50	196.3			100								1.15	250	
25		HiScale 10	15.7	20% ethanol with 0.4 M NaCl	45–55	16	300	N/A	N/A	1.0	N/A	14	1070			
		XK 16 and HiScale 16	40.2			10								1.1	15	
		XK 26 and HiScale 26	106.1			27								1.13	40	
35		XK 50 and HiScale 50	392.5	20% ethanol with 0.4 M NaCl	45–55	10	300	N/A	N/A	1.1	N/A	150	450			
		HiScale 10	19.6			14								1.0	14	1070
		XK 16 and HiScale 16	70.3			10								1.06	8.6	
35	XK 26 and HiScale 26	185.7	20% ethanol with 0.4 M NaCl	45–55	27	300	N/A	N/A	1.1	N/A	23	260				
	XK 50 and HiScale 50	686.9			100								1.06	86		

¹ Part of the CDM program from GE. Packing protocols for CDM-developed resins are not verified but are based on packing protocols for standard products having similar pressure-flow characteristics.

² The PF values in this table are specific for the lab-packing method used. The PF shown here bears no relation to the PF described for large-scale column packing procedures.

³ d_{50v} = median particle size of the cumulative volume distribution.

Capto Core 700 and Capto MMC

Packing parameters

Your parameters				Resin preparation		Column packing						Post-column packing								
Resin	Packed bed height (cm)	Column type	Bed volume (mL)	Slurry/packing solution	Slurry concentration (%)	Settling flow The flow that consolidates the resin suspension (slurry)		Packing flow The flow that compresses the bed		Compression Bed is further compressed by turning the adapter manually into the bed		Conditioning flow The flow to get a uniform bed over the entire column height								
						mL/min	cm/h	mL/min	cm/h	Mechanical compression packing factor (PF) ¹	Flow compression (mm)	mL/min	cm/h							
Capto Core 700 Particle size d_{50v}^2 : 85 μ m	10	Tricorn 5	2.0	20% ethanol with 0.2 M NaCl	45-55	2.29	700	2.29	700	N/A	1	N/A	N/A							
		Tricorn 10	7.9			9.16		9.16												
		XK 16 and HiScale 16	20.1			10		N/A						N/A	1.15	N/A	25			
		XK 26 and HiScale 26	53.1			27									1.15		66			
		XK 50 and HiScale 50	196.3			100									1.15		250			
		XK 16 and HiScale 16	40.2			10									1.1		15			
	XK 26 and HiScale 26	106.1	27	1.1	40															
	XK 50 and HiScale 50	392.5	100	1.1	150															
	20	35	XK 16 and HiScale 16	70.3	20% ethanol with 0.4 M NaCl	45-55	10	300	N/A	N/A	1.06	N/A	8.6	260						
			XK 26 and HiScale 26	185.7			27				1.06		23							
			XK 50 and HiScale 50	686.9			100				1.06		86							
		10	Tricorn 5	2.0			10 mM NaCl				40-60		1.8		540	9.8	3000	N/A	1	N/A
Tricorn 10			7.9	7.1									39.3							
XK 16 and HiScale 16			20.1	25									N/A			N/A				
XK 26 and HiScale 26	53.1	66	1.15	66																
XK 50 and HiScale 50	196.3	250	1.15	250																
20	XK 16 and HiScale 16	40.2	25	1.1	25															
	XK 26 and HiScale 26	106.1	66	1.1	66															
	XK 50 and HiScale 50	392.5	250	1.15	250															
35	XK 16 and HiScale 16	70.3	25	1.02	14															
	XK 26 and HiScale 26	185.7	66	1.03	37															
	XK 50 and HiScale 50	686.9	250	1.03	140															

¹ The PF values in this table are specific for the lab-packing method used. The PF shown here bears no relation to the PF described for large-scale column packing procedures.

² d_{50v} = median particle size of the cumulative volume distribution.

Capto adhere and Capto adhere ImpRes

Packing parameters

Your parameters				Resin preparation		Column packing						Post-column packing							
Resin	Packed bed height (cm)	Column type	Bed volume (mL)	Slurry/packing solution	Slurry concentration (%)	Settling flow The flow that consolidates the resin suspension (slurry)		Packing flow The flow that compresses the bed		Compression Bed is further compressed by turning the adapter manually into the bed		Conditioning flow The flow to get a uniform bed over the entire column height							
						mL/min	cm/h	mL/min	cm/h	Mechanical compression packing factor (PF) ¹	Flow compression (mm)	mL/min	cm/h						
Capto adhere Particle size d_{50v}^2 : 75 μm	10	Tricorn 5	2.0	10 mM NaCl	40-60	0.25	76	10	3000	N/A	1	N/A	N/A						
		Tricorn 10	7.9			1		40											
		HiScale 10	7.9	20% ethanol with 0.4 M NaCl	60	20	1528	20	1528	1.0	N/A	20	1528						
		XK 16 and HiScale 16	20.1	20% ethanol with 0.4 M NaCl	45-55	25	750	N/A	N/A	1.1	N/A	25	750						
		XK 26 and HiScale 26	53.1			66				1.15		66							
		XK 50 and HiScale 50	196.3			250				1.15		250							
	20	HiScale 10	15.7	20% ethanol with 0.4 M NaCl	45-55	20	1528	20	1528	1.0	N/A	20	1528						
		XK 16 and HiScale 16	40.2			25				750		N/A		N/A	1.1	25			
		XK 26 and HiScale 26	106.1			66									1.1	66			
		XK 50 and HiScale 50	392.5			250									1.15	250			
		25	HiScale 10			19.6				60		20		1528	20	1528	1.0	20	1528
		35	XK 16 and HiScale 16			70.3				45-55		750		N/A	N/A	1.02	14		
	XK 26 and HiScale 26		185.7	66	1.03	37													
	XK 50 and HiScale 50		686.9	250	1.03	140													
	Capto adhere ImpRes Particle size d_{50v}^2 : 40 μm	10	Tricorn 5	2.0	10 mM NaCl	45-55	7.4	2250	7.4	2250	N/A	1	N/A	N/A					
Tricorn 10			7.9	29.4			29.4												
HiScale 10			7.9	20% ethanol with 0.4 M NaCl	45-55	10.3	783	14.5	1108	1.05	N/A	14.5	1108						
XK 16 and HiScale 16			20.1			20				600		N/A		N/A	1.12	20			
XK 26 and HiScale 26			53.1			53									1.12	53			
XK 50 and HiScale 50			196.3	196	1.12	196													
20		HiScale 10	15.7	20% ethanol with 0.4 M NaCl	45-55	10.3	783	14.5	1108	1.05	N/A	14.5	1108						
		XK 16 and HiScale 16	40.2			17				500		N/A		N/A	1.12	17			
		XK 26 and HiScale 26	106.1			44									1.12	44			
		XK 50 and HiScale 50	392.5			164									1.12	164			
		25	HiScale 10			19.6				63		10.3		783	14.5	1108	1.05	14.5	1108
		35	XK 16 and HiScale 16			70.3				45-55		350		N/A	N/A	1.09	12		
XK 26 and HiScale 26			185.7	31	1.09	31													
XK 50 and HiScale 50			686.9	115	1.09	115													

¹ The PF values in this table are specific for the lab-packing method used. The PF shown here bears no relation to the PF described for large-scale column packing procedures.

² d_{50v} = median particle size of the cumulative volume distribution.

Capto MMC ImpRes, Capto Phenyl ImpRes, and Capto Butyl ImpRes

Packing parameters

Your parameters				Resin preparation		Column packing						Post-column packing							
Resin	Packed bed height (cm)	Column type	Bed volume (mL)	Slurry/packing solution	Slurry concentration (%)	Settling flow The flow that consolidates the resin suspension (slurry)		Packing flow The flow that compresses the bed		Compression Bed is further compressed by turning the adapter manually into the bed		Conditioning flow The flow to get a uniform bed over the entire column height							
						mL/min	cm/h	mL/min	cm/h	Mechanical compression packing factor (PF) ²	Flow compression (mm)	mL/min	cm/h						
Capto MMC ImpRes Particle size d_{50v} ³ : 40 μ m	10	Tricorn 5	2.0	10 mM NaCl	45-55	7.5	2250	7.5	2250	N/A	1	N/A	N/A						
		Tricorn 10	7.9			29.4		29.4											
		XK 16 and HiScale 16	20.1			20		600						1.15	N/A	20	600		
	XK 26 and HiScale 26	53.1	53	1.15	53														
	XK 50 and HiScale 50	196.3	196	1.15	196														
	20	XK 16 and HiScale 16	40.2	20% ethanol with 0.4 M NaCl	35-45	17	500	N/A	N/A	1.12	N/A	17	500						
		XK 26 and HiScale 26	106.1			44								1.12	44				
		XK 50 and HiScale 50	392.5			164								1.12	164				
	35	XK 16 and HiScale 16	70.3	20% ethanol with 0.4 M NaCl	35-45	12	350	N/A	N/A	1.09	N/A	12	350						
		XK 26 and HiScale 26	185.7			31								1.09	31				
		XK 50 and HiScale 50	686.9			115								1.09	115				
	Capto Phenyl ImpRes ¹ Capto Butyl ImpRes ¹ Particle size d_{50v} : 40 μ m	10	Tricorn 5	2.0	20% ethanol	45-55	3.8	1125	3.8	1125	N/A	1	N/A	N/A					
Tricorn 10			7.9	14.7			14.7												
HiScale 10			7.9	10.3			783		16						1222	1.05	N/A	16	1222
20		XK 16 and HiScale 16	20.1	20% ethanol	45-55	20	600	N/A	N/A	1.15	N/A	20	600						
		XK 26 and HiScale 26	53.1			53								1.15	53				
		XK 50 and HiScale 50	196.3			196								1.15	196				
25		HiScale 10	15.7	20% ethanol with 0.4 M NaCl	63	10.3	783	16	1222	1.05	N/A	16	1222						
		XK 16 and HiScale 16	40.2			17								1.12	17				
		XK 26 and HiScale 26	106.1			44								1.12	44				
35		XK 50 and HiScale 50	392.5	20% ethanol	45-55	164	500	N/A	N/A	1.12	N/A	164	500						
		HiScale 10	19.6			10.3								783	14	1070	1.05	14	1070
		XK 16 and HiScale 16	70.3			12								350	N/A	N/A	1.09	12	350
XK 26 and HiScale 26	185.7	31	1.09	31															
XK 50 and HiScale 50	686.9	115	1.09	115															

¹ Part of the CDM program from GE. Packing protocols for CDM-developed resins are not verified but are based on packing protocols for standard products having similar pressure-flow characteristics.

² The PF values in this table are specific for the lab-packing method used. The PF shown here bears no relation to the PF described for large-scale column packing procedures.

³ d_{50v} = median particle size of the cumulative volume distribution.

Capto Phenyl (High Sub), Capto Butyl, Capto Octyl, KappaSelect, LambdaFabSelect, IgSelect, Alpha-1-Antitrypsin Select, Capto Chelating, Capto Blue, Capto Blue (High Sub), CaptoDeVirs, and VII Select

Packing parameters

Your parameters				Resin preparation		Column packing						Post-column packing	
Resin	Packed bed height (cm)	Column type	Bed volume (mL)	Slurry/packing solution	Slurry concentration (%)	Settling flow The flow that consolidates the resin suspension (slurry)		Packing flow The flow that compresses the bed		Compression Bed is further compressed by turning the adapter manually into the bed		Conditioning flow The flow to get a uniform bed over the entire column height	
						mL/min	cm/h	mL/min	cm/h	Mechanical compression packing factor (PF) ²	Flow compression (mm)	mL/min	cm/h
Capto Phenyl (High Sub)	10	Tricorn 5	2.0	20% ethanol with 0.2 M NaCl	45-55	1	300			1.05		2	600
Capto Butyl		Tricorn 10	7.9			4				1.08		8	
Capto Octyl ¹		XK 16 and HiScale 16	20.1			25				1.1		25	
KappaSelect ¹		XK 26 and HiScale 26	53.1			66				1.15		66	
Lambda-FabSelect ¹		XK 50 and HiScale 50	196.3			250				1.15		250	
Alpha-1-Antitrypsin Select ¹	20	XK 16 and HiScale 16	40.2	20% ethanol with 0.4 M NaCl	45-55	25	750	N/A	N/A	1.1	N/A	25	750
Capto Chelating ¹		XK 26 and HiScale 26	106.1			66				1.1		66	
Capto Blue		XK 50 and HiScale 50	392.5			250				1.15		250	
Capto Blue (High Sub) ¹	35	XK 16 and HiScale 16	70.3			25				1.02		14	420
CaptoDeVirs ¹		XK 26 and HiScale 26	185.7			66				1.03		37	
VII Select ¹		XK 50 and HiScale 50	686.9			250				1.03		140	

¹ Part of the CDM program from GE. Packing protocols for CDM-developed resins are not verified but are based on packing protocols for standard products having similar pressure-flow characteristics.

² The PF values in this table are specific for the lab-packing method used. The PF shown here bears no relation to the PF described for large-scale column packing procedures.

³ d_{50v} = median particle size of the cumulative volume distribution.



gelifesciences.com

GE, the GE Monogram, Axichrom, Capto, Chromaflow, HiScale, MabSelect, MabSelect SuRe, and Tricorn are trademarks of General Electric Company.
© 2019 General Electric Company.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them.
A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.

GE Healthcare Bio-Sciences AB, Björkgatan 30, 751 84 Uppsala, Sweden

GE Healthcare Bio-Sciences Corp., 100 Results Way, Marlborough, MA 01752, USA

GE Healthcare Europe GmbH, Munzinger Strasse 5, D-79111 Freiburg, Germany

GE Healthcare Japan Corp., Sanken Bldg., 3-25-1, Hyakunincho Shinjuku-ku, Tokyo 169-0073, Japan

GE Healthcare UK Ltd., Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA, UK

For local office contact information, visit gelifesciences.com/contact

KA4551032119SG